



USER GUIDE

Apollo[®] 20 mL PES & PP .45 μ m High-Performance Centrifugal Microfilters

Note: This product is offered for research use only. Not for clinical use, diagnostic procedures, or for preparation of fluids to be used for human injection. This is not a sterilizing filter. This is a pilot production product that is in final development. Preliminary applications data are available, and a Quality Control System is in place for this product.

Apollo 20 mL MF PES and PP concentrators are disposable single-use microfiltration devices for the removal, concentration or purification of suspended particulates and colloids. They are far superior to alternatives in combined simplicity, speed, capacity and recovery. This is due to their unique conical design (US Patent 6,269,957, US Patent 6,357,601, PCT patents pending), providing a high ratio of membrane area to sample size. This, in turn, provides a high degree of concentration in a single spin as well as better control of polarization and fouling at the membrane surface. Apollo has the largest available sample volumes for a given centrifuge tube size.

SPECIFICATIONS Volumes

	<u>Maximum Sample</u>
34° angle rotor:	13.5 mL
With swing-head rotor:	20 mL

Total Volume In Concentrator & Filtrate Tube

<u>Swing Head</u>	<u>34° Angle</u>	<u>Resulting Deadstop</u>
NA	<= 14 mL	13 μ L* / 28 μ L**
<= 20.0 mL	20.6 mL	67 μ L
20.6 mL	21.3 mL	100 μ L
22.0 mL	22.4 mL	200 μ L
24.7 mL	24.4 mL	500 μ L
27.2 mL	26.2 mL	1 mL
28.9 mL	27.4 mL	1.5 mL
30.2 mL	28.5 mL	2 mL

* With any port oriented outboard

** With any port oriented inboard

Dimensions

Active membrane area=12 cm²
Collection tube: OD=29.2 mm,
Length (incl. cap)=118.2 mm
Filter: Length (filter tip to top flange)=
68 mm OD (below top flange)=26.8
mm

Maximum Centrifugal Force

With swing-head rotor: 500 to 2000 rcf
typical rotor speed, not to exceed 4500 rcf

Materials

Membrane: PES=Hydrophilic
polyethersulfone, PP=polypropylene
Concentrator, collection tube and cap:
Polypropylene

Environmental Resistance

Temperature: 34.7 °C, 120 °F, max.
Do not autoclave
Limit of pH: 1 to 14

HOW TO USE THIS PRODUCT

Preparations

Make sure it will fit in your centrifuge

Prepare a 50 mL carrier accepting a 118 mm length tube in centrifuge. Either fixed angle or swing head rotors can be used, although performance is better in a swing head. Check clearance of tube to both swing mechanism and rotor cover or centrifuge lid.

Suggested protocol if removal of trace metals is required

Rinse inside surface with a squirt bottle of 1% nitric acid and pour out. Fill device with ~20 mL 1% nitric acid and spin for 2 to 3 minutes to flush down to a retentate volume <0.1 mL (see figure). Pour out dilute acid from tube and filter. Rinse and empty tube and filter device 2 to 3 times with deionized water. Shake out water droplets and add sample.

Operation

1. Add sample and cap tube snugly.

An internal vent hole near the lip permits air from the collection tube to pass into the concentrator to maintain maximal flow without release of aerosols. *Avoid scraping membrane when filling.*

2. Place assembly into rotor. Counterbalance with a similar device or tube of the same weight and spin. *Note specified centrifugal force limits and observe maximum relative centrifugal force rating for the rotor.*

3. Spin for the required time (see typical performance graph for PES)

Spin at the suggested speed to achieve the desired concentration factor. Exceeding the maximum centrifugal force limits, specified above, may cause retentate leakage. *Avoid excessive rcf, improved selectivity maybe achieved by avoiding excessive rcf.*

4. Harvest retentate or filtrate

For best recovery, remove retentate in <10 min. Upon standing, wicking by the spun, partly desiccated membrane can cause continued filtration, further reducing retentate volume. For retentate volumes <50 μ L, mass recovery is improved by adjusting volume with buffer to about 50 μ L before recovery, and/or by subsequently adding 50-200 μ L of buffer to the device, mixing by vortexing or aspirating into and out of the tip several times and recovering the wash as well.

Chemical Compatibility for Common chemicals (✓ = acceptable; **X** = *not recommended*, ? = insufficient data)

Acids and Bases


	PES	PP		PES	PP		PES	PP
Acetic acid 10%	✓	✓	Hydrochloric acid 35%	✓	✓	Ammonium hydroxide 11.4%	✓	✓
Acetic acid 30%	✓	✓	Nitric acid 27	?	✓	Potassium hydroxide 15%	✓	✓
Acetic acid 90%	✓	✓	Nitric acid 67%	X	✓	Sodium hydroxide 11%	✓	✓
Glacial acetic acid	✓	✓	Sulfuric acid 16%	?	✓	Sodium hydroxide 22%	✓	✓
Hydrochloric acid 3.3%	✓	✓	Sulfuric acid 96%	X	✓			
Hydrochloric acid 20%	✓	✓	Ammonium hydroxide 5.7%	✓	✓			

Organic Solvents, Miscellaneous Chemicals

	PES	PP		PES	PP		PES	PP
Acetone	X	✓	Dimethyl sulfoxide	X	✓	Methanol	✓	✓
Acetonitrile	✓	✓	Ethanol	✓	✓	Methyl acetate	X	✓
Amyl acetate	✓	✓	Ethyl acetate	X	✓	Methylene chloride	X	✓
Amyl alcohol	X	✓	Ethyl ether	✓	✓	Methyl ethyl ketone	X	✓
Benzene	✓	✓	Ethylene dichloride	X	✓	Methyl isobutyl ketone	✓	✓
Benzyl alcohol	X	✓	Ethylene glycol	✓	✓	Peanut oil	✓	✓
Butanol	✓	✓	Formaldehyde 4%	✓	✓	Propylene glycol	✓	✓
Butyl acetate	✓	✓	Formaldehyde 37%	✓	✓	Pyridine	X	✓
Carbon tetrachloride	✓	✓	Glycerol	✓	✓	Tetrachloroethylene	✓	✓
Cellosolve acetate	✓	✓	Hexane, dry	X	✓	Toluene	✓	✓
Chloroform	X	✓	Isopropanol	✓	✓	Tetrahydrofuran	X	✓
Cotton seed oil	✓	✓	Isopropyl acetate	✓	✓	Tetrahydrofuran/water 50%v	?	✓
Cyclohexanone	X	✓	Kerosene	✓	✓	Xylene	✓	✓
Dimethyl formamide	X	✓						

Some of the recommended chemicals listed above may affect membrane performance, thereby altering the recoveries, passage, and /or spin times.

TYPICAL PERFORMANCE (Actual conditions will vary with details of initial solution temperature, concentration, and sample characteristics.)



Note: Apollo 20 PP .45µm is a hydrophobic membrane and will normally require a minimum 2000 rcf to overcome the breakthrough pressure and develop a reasonable flow for most aqueous samples. (spun 5 minutes @ 2000rcf with MilliQ Biocel water = ~ 18 mL of filtrate) Reasonable flows can be achieved using lower rcf if the membrane is pre-wet with an organic solvent such as isopropanol.

Technical Assistance

Call, fax, or e-mail us at the numbers below for help. You can also visit us on the Internet www.orbio.com for the most up-to-date technical information on the Apollo family of products.

To reorder, please refer to Product Numbers below.

Product Name	Rating	Identification	Qty/Pk	Order No.
Apollo 20 mL PES .45µm	.45µm	Sample pack, sealed devices and User Guide	2 ea.	AP20PS4500
Apollo 20 mL PES .45µm	.45µm	Rack of filters in capped 50 mL tubes	25 ea.	AP20PS4510
Apollo 20 mL PES .45µm	.45µm	10 racks, bulk bags of filters, 50 mL tubes, caps	250 ea.	AP20PS4520
Apollo 20 mL PP .45µm	.45µm	Sample pack, sealed devices and User Guide	2 ea.	AP20PP4500
Apollo 20 mL PP .45µm	.45µm	Rack of filters in capped 50 mL tubes	25 ea.	AP20PP4510
Apollo 20 mL PP .45µm	.45µm	Rack of filters in capped 50 mL tubes	250 ea.	AP20PP4520
50 mL tubes and cap	-	Rack of 25 ea tubes and caps for Apollo 20 mL	1 ea.	AP2000000
50 mL tubes and cap	-	Case of tubes & caps for Apollo 20 mL	500 ea.	3191-870-00

Orbital Biosciences, LLC

4 Winsor Lane • P.O. Box 304 • Topsfield, MA 01983 • USA
Tel. 978-887-5077 • FAX 978-887-2906 • e-mail: orbio@shore.net

Since Orbital Biosciences has no control over the actual use conditions of its products, it shall not be liable for any consequential damages or economic loss or property damage incurred by Buyer. In connection with such use, Buyer's remedies under Orbital Bioscience's warranty shall be limited to replacement of the product or component, which failed to perform according to specifications herein stated. AP20 PES & PP USER GUIDE, REV 1, 04/22/02